



QUICK START GUIDE

Welcome to the Pegasus 3D Footscan experience

POWERED BY TECHMED  3D



Overview – True 3D iPad Scanning

The Pegasus 3D Footscan application (app) is an integrated app used to capture true 3d footscan anatomy of the patient's feet.

The APP facilitates instant, digital transmission of the footscans along with an interactive, digital order form to the Pegasus lab for orthotic design and fabrication.

A confirmation email is automatically generated from the lab to confirm successful receipt of footscans and order form.

Prerequisites

To use the Pegasus 3D Footscan app, the following items are needed:

- _ iOS Device (iPad or iPhone)
- _ Structure Sensor™ scanner
- _ Internet Access and an e-mail account on the iOS device so that the captured scan files and completed order forms can be e-mailed to Pegasus
- _ The Pegasus 3D Footscan app (available on the App Store).



Installation

Visit the App Store and download the Pegasus 3D Footscan app.

[Click here to learn how to install your app >](#)

[Account setup and subscription >](#)

Scanning Process

[Click here to learn how to scan a plantar surface >](#)

Sending Orders

E-mail orders to:
3dfootscan@pegasusorthotics.com

Other Resources

[Backup procedures >](#)

[Delete files or patients >](#)

[Precision Assessment procedure >](#)

Support

Pegasus Orthotics Ltd
support@pegasusorthotics.com
+44 (0)28 9210 7200 ext. 201

